

# ZEMITS SKIN ANALYSIS

For professional skin assessment


Comprehensive skin  
analysis across 12 skin  
concerns



# ZEMITS SKIN ANALYSIS

 New prospects for your aesthetic practice

The Zemits skin analyser is an advanced device specially designed for detailed skin analysis before, during, and after skin treatments. It utilises artificial intelligence and RGB, UV, and PL technologies to provide fast and accurate results in just 15 seconds



Large screen 10.1



Accurate and quick analysis of results



Futuristic design



Expansion  
16 megapixels



# ADVANTAGES

★ **Modern Technologies:** combines RGB, UV and PL technologies with artificial intelligence, creating 3D images for deep skin analysis.


★ **Accurate diagnosis:** A Professional Analyser provides detailed information about the skin's condition, allowing you to identify problem areas and choose optimal care.

★ **Individual approach:** Each client receives a unique care plan that takes into account their specific needs and skin characteristics.


★ **Proactive Care:** Regular skin analysis is not just about monitoring changes; it's about being in control. It allows you to adjust care programs to achieve the best results, making you feel proactive and in charge of your skin's health.

## COMPREHENSIVE SKIN DIAGNOSTIC

✓ ACROSS 12 CONCERN WITH ONE SCAN



Instant **skin analysis**  
powered by AI!



► Zemits Skin  
► Analysis

10.1

Large screen for  
ultra-precise analysis:  
10.1 inches

16

Screen expansion of 16  
megapixels - will not miss a  
single wrinkle

15

Just 15 seconds to analyse and  
forecast the client's skin condition  
up to 3-5 years ahead


12

Comprehensive analysis across 12  
skin factors - the degree of skin  
hydration, the presence of acne,  
the state of the skin texture and  
enlarged pores

# OPEN NEW SOURCES OF INCOME



with Skin Analyser



## Up-sale of services and skincare:

The skin analyser enables you to develop personalised skin care plans that include suggestions for specific products and treatments the aesthetician provides. This creates opportunities to upsell extra services and retail products that cater to the client's requirements. With personalised recommendations, practitioners can propose additional treatments, like facial massages and specialised masks, improving skin treatment results and boosting the salon's profitability.

Prompt skin analysis and treatment plan.



## PORES

RGB spectral visualization of pores: analyzes pores as low as 0.02-0.05 mm in diameter.



## WRINKLES

RGB Spectral visualization of wrinkles: analysis of the depth of skin creases and its irregularities. The condition of wrinkles is accurately determined by the detailed analysis.



## ACNE

RGB spectral visualization of acne: analysis of dark spots on round-shaped skin and darker than normal colors.



## UV PIGMENTATION

Pigmentation analysis via UV spectral imaging: This is a display not only of the current pigmentation situation, but also of the prediction of future formations. The darker areas require more attention.




## TEXTURE

Spectral visualization of the PL skin texture: This is the display of the relief, texture of the skin surface, as well as the prediction of future age-related changes.




## UV PIGMENTATION

UV damage: Pigmentation can occur at different levels of the skin. It is formed under the negative action of ultraviolet rays and is a stain. The device fixes the depth and size of age spots.



## SENSITIVE AREA

Sensitive areas: Fixation of the zonal sensitivity of the skin, its susceptibility to allergic reactions. This can be both a seasonal allergy and a reaction to cosmetics with harmful elements in the composition



## UV PORPHYRIN

Spectral visualization of porphyrin UV: Sebum secretion and the area of distribution of black dots.



## UV PORPHYRIN

Spectral visualization of porphyrin UV: Sebum secretion and the area of distribution of black dots.



## UV DAMAGE

UV damage: Pigmentation can occur at different levels of the skin. It is formed under the negative action of ultraviolet rays and is a stain. The device fixes the depth and size of age spots.



## UV MOISTURE

UV hydration: How much the moisture level of the dermis corresponds to the norm.



## SKIN AGING

Skin aging: Shows the dynamics of the aging process of the skin, provided that the client will not change anything in the care.

# TESTIMONIALS



**thrivekalamazoo**

Thrive Health & Wellness

Zemits skin analyzer provides a detailed analysis of your skin, giving valuable insights about its condition and needs. With this information, our expert team can recommend a personalized skin care routine tailored specifically for clients! 😊😊



[READ MORE](#)



**FBE**

NYC Master Aesthetician

This advanced system is designed to analyse facial skin, capturing, displaying, and comparing results. It can identify a wide range of skin issues, including age spots, pores, moisture levels, texture, and wrinkles. With this comprehensive analysis, we can provide personalised recommendations for skincare products, treatments, and more, all aimed at giving your skin a radiant glow!



[READ MORE](#)

# TESTIMONIALS



Emma Lee

Skin & Aesthetics

Are you struggling with wrinkles, fine lines, ageing skin or saggy skin?

Let me help you get your skin looking and feeling great!  
With Zemits Skin Analysis I will provide a skin analysis consultation, advice & a taster facial



[READ MORE](#)



VZ Beautytherapist

ESHER&EPSOM ❤

Skin Analyser is an advanced imaging system that uses a high resolution digital camera with 3 spectrum imaging to analyse 10 facial problems including acne, open and blocked pores, hyper-pigmentation, broken capillaries, uv damage and wrinkles. It produces detailed images and reports.



[READ MORE](#)

# TECHNICAL SPECIFICATION

**Power:** 45 W

**High-voltage:** 110-220 V

**Pixel:** 16 megapixels

**Screen size:** 10.1 inches

**Current:** 0.2 A, 50 Hz

**Device size:** 42x36x56 cm

**Housing size:** 50x45x67 cm

**Weight:** 8 KG / 12 KG



# Why Zemits?

## Financing

Acquire our cutting-edge device through our in-house instalment program. Start benefiting today and defer your payments until a later date.

## Marketing support

With every Zemits device we provide comprehensive marketing support to help expand your clientele.

## Complementary Professional training with every device

In-depth theory and practical training delivered face-to-face or virtually is included with every device purchase.

## 24-month warranty and UK service centre

Rest assured of the safety and reliability of your device.



**Join Zemits Family**

Discounts   Offers   Gifts

# VIRTUAL DEMO

Guarantee of confidence  
in the purchase

## What is a VIRTUAL DEMO?

A new, unique, absolutely free, but very useful service.

Go through an overview of the selected device, including treatment options and protocols and have face-to-face consultations with Zemits certified Aesthetic Sales advisers to get information about the best technology, consider clientele demand and your experience for the best investment with high ROI.



## VIRTUAL DEMO

Book Demo



Click to book

# OUR CONTACTS

We invite you to our showroom in Cheam, SW of London, to ensure the quality of your device and meet our fantastic team of experts! Booking in advance is required.

 020 4532 5922

## Our Business Hours:

09:00 - 18:00 (GMT) MON-FRI,

Place an Order Online 24/7